

MYCO PNEUMO ONE®

System for the presumptive identification and antibiotic susceptibility test of *Mycoplasma pneumoniae*

1. INTRODUCTION

Mycoplasma pneumoniae causes up to 10-40% of community-acquired pneumonia. The incidence of *M. pneumoniae* pneumonia is higher among children and young adults.

The symptoms of *M. pneumoniae* infections of the upper and lower respiratory tract are generally mild and often self-limiting. The most frequent extrapulmonary complications are associated with the CNS, the heart and the skin. Skin diseases usually occur with erythematous or transient maculopapular blisters, but can sometimes develop into Stevens-Johnson syndrome.

A percentage of asthmatic exacerbations is related to *M. pneumoniae*, as a consequence of chronic diseases of the lower respiratory tract. Numerous studies have also shown a link between *M. pneumoniae* and recurrent asthmatic processes or that do not respond to conventional treatments.

The presence of these microorganisms has also been associated with ocular infections⁽¹⁻⁸⁾.

As a result, the logical procedure would be to eliminate these pathogens as soon as possible using antibiotics. However, antibiotics are not recommended for the treatment of acute asthma exacerbations, except when necessary. These differences underscore the need to define the role of antimicrobials that are active against atypical primarily macrolide pathogens, but also tetracyclines and fluoroquinolones in the treatment of asthma⁽⁵⁻⁷⁾.

A system that allows accelerated growth, with antimicrobial susceptibility in just 48 hours, without additional equipment, can be a useful tool in the hands of the microbiologist and the clinician in the diagnosis of *Mycoplasma pneumoniae*.

2. PRINCIPLE

System consisting of a plate containing 16 conical wells designed for better visualization of colorimetric reactions that occurs as a result of growth and antimicrobial susceptibility testing in culture media especially formulated for selective culture and AST of *Mycoplasma pneumoniae*.

The plate also contains control wells for the detection of cross-reactions and / or non-specific reactions caused by the type of sample, inadequate manipulation, contamination.

The system is inoculated with the test sample and incubated at 36 ± 1 ° C for 48 hours. Following incubation, the micro-organism (or the eventual copresence of several micro-organisms) is identified by combining the colorimetric reactions that take place in the wells.

Antimicrobial susceptibility is designed to help in the identification and evaluation of cut off susceptibility and resistance of antibiotics as recommended by CLSI.

Microbiological diagnoses can be confirmed by serological tests, microscopic observation, culture or molecular techniques directly from the positive wells.

3. CONTENT OF MYCO PNEUMO ONE KIT REF. MS01352

10 Identification panels

10 Storage Liquide Solutions *Mycoplasma spp.*

4. REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

Culture Medium Plates enriched for *Mycoplasma spp.*


Reagents for Gram Stain

Sterile liquid paraffin

Sterile swabs

Generic laboratory material

5. PANEL IDENTIFICATION DESCRIPTION



WELL	DESCRIPTION	CONTENT
1	<i>M.pneumoniae/M.spp</i>	Selective medium for the presumptive identification of <i>Mycoplasma spp.</i>
2	<i>M.pneumoniae/M.spp</i>	Selective medium for the presumptive identification of <i>Mycoplasma spp.</i>
3	<i>M.pneumoniae/M.spp</i>	Selective medium for the presumptive identification of <i>Mycoplasma spp.</i>
4	<i>M.pneumoniae/M.spp</i>	Selective medium for the presumptive identification of <i>Mycoplasma spp.</i>
5	<i>M.pneumoniae/M.spp</i>	Selective medium for the presumptive identification of <i>Mycoplasma spp.</i>
6	Cross Reaction	Selective medium for the determination of cross-reactions.
7	Fungal Growth	Selective medium for the identification of fungi and yeasts in the sample
8	Bact Growth	Selective medium for the identification of bacteria not belonging to the <i>Mycoplasma</i> species.
9	LEV 1	Selective medium for <i>Mycoplasma spp.</i> containing Levofloxacin 1 µg/mL
10	MXF 0,5	Selective medium for <i>Mycoplasma spp.</i> containing Moxifloxacin 0,5 µg/mL
11	TET 2	Selective medium for <i>Mycoplasma spp.</i> containing Tetraciclina 2 µg/mL
12	E 0,5	Selective medium for <i>Mycoplasma spp.</i> containing Erythromycin 0,5 µg/mL
13	E 1	Selective medium for <i>Mycoplasma spp.</i> containing Erythromycin 1 µg/mL
14	AZ 0,5	Selective medium for <i>Mycoplasma spp.</i> containing Azitromicina 0,5 µg/mL
15	AZ 1	Selective medium for <i>Mycoplasma spp.</i> containing Azitromicina 1 µg/mL
16	AZ 16	Selective medium for <i>Mycoplasma spp.</i> containing Azitromicina 16 µg/mL

6. EXECUTION OF THE TEST

6.1 SAMPLES AND PROCEDURE

SAMPLES:

Pharyngeal swab, nasopharyngeal swab, bronchial secretions, tracheal aspiration, pleural fluid aspiration, transthoracic needle aspiration, trans-tracheal needle aspiration, broncho-alveolar lavage, pulmonary needle aspiration, abscess needle aspiration, blood cultures.

To obtain the best performances with the use of this system, the sample must be collected aseptically according to the methodology implemented in each hospital center and before starting the antibiotic treatment.

The sample must be tested immediately after collection (Paragraph 6.3).

Do not refrigerate or freeze. Protect from light and sun.

The **MYCO PNEUMO ONE**® system must be used by a qualified microbiologist.

In case of doubt in the procedure or in the interpretations, contact the product specialist of C.P.M. on site or directly at the central office.

Whenever possible, it is recommended that respiratory tract samples are grown in conventional media such as Columbia Agar with 5% sheep blood.

6.2 PROCEDURE FOR TEST EXECUTION

In the case of pharyngeal and nasopharyngeal exudates, resuspend the swab containing the sample in the vial of the Storage Liquide Solutions *Mycoplasma spp.* provided in the kit, shake the swab and press it against the walls to obtain a homogeneous suspension. In the case of other types of respiratory samples to be analyzed, using a Pasteur pipette or a sterile swab, resuspend a portion of the sample directly in the Storage Liquide Solutions *Mycoplasma spp.*

Of the suspension obtained, inoculate 150 µL in each well of the **MYCO PNEUMO ONE** panel.

WARNING the samples must be collected according to the sample collection procedure established in each laboratory and based on the standard guidelines ⁽¹⁷⁻²⁰⁾ (see section 6.3).

/! Add two drops of sterile paraffin to wells 1 to 16.

/! The other samples of the tract respiratory to be analyzed using the **MYCO PNEUMO ONE** system must be processed as established in each laboratory, according to what is reported in section 6.3.

The microbiologist must evaluate what type of sample should be inoculated directly into the Storage Liquide Solution *Mycoplasma spp.*

Incubation

Incubate at 36 ± 1 ° C for 48 hours, in a traditional laboratory incubator.

The incubation can be extended up to 5 days, taking into account the time required for the growth of the microorganisms included in the identification panel and their concentration in the sample in examination.

6.3 RECOMMENDATIONS FOR THE COLLECTION AND SAMPLE ANALYSIS

Each hospital and / or laboratory center must follow the procedures established internally for the selection of the samples to be studied. The recommendations presented here are for guidance only .

Sample collection, aspects to consider: ^{(1) (4) (5)}

- Assess the risk / benefit of patient sample collection, especially in case of invasive samples.
- Samples must be transported in suitable containers.
- Sample collection must be carried out under aseptic conditions
- An adequate amount of sample must be collected.

if it was necessary to transport or store the sample until its analysis, it is recommended to use an adequate transport medium. Use of transport media such as Amies or Stuart interfere with the results of this kit.

The sample must be analyzed before two hours from collection, and conserved into the Storage Liquide Solution *Mycoplasma spp.* provided in the kit, which contains a suitable composition for preserve the conservation of fastidious microorganisms such as those detectable with this kit.

The samples of the respiratory tract should not be stored in conventional media for more than 24 hours ⁽¹⁸⁻²⁰⁾ , and in these conditions are **NOT** recommended for use with this kit.

Sample processing: The sample must be handled according to the procedure established in each laboratory and must be prepared, stained and inoculated properly taking in account the type of sample received, its clinical diagnosis and the physician's prescription. Depending on the type of sample, the laboratory should evaluate and determine the need for centrifugation, homogenization or other procedure before its inoculation.

The Storage Liquide Solution *Mycoplasma spp.* allows the conservation of the sample at room temperature for 24 hours, 3 days at 4 ° C, or freezing at -20 / -70 ° C ⁽¹⁸⁻²⁰⁾. Although the samples of the respiratory tract can sometimes be polymicrobial, the Storage Liquide Solution *Mycoplasma spp.* contains inhibitors and is selective for *Mycoplasma spp.*, It is up to the microbiologist to use the advantages of the Storage Liquide Solution *Mycoplasma spp.* for further studies.

It is recommended to analyze the results obtained with MYCO PNEUMO ONE together with the traditional methods and complementary tests established in each laboratory. This kit combines the utility of the use of selective and / or chromogenic media with the possibility of carrying out any type of confirmation using the contents of the well where bacterial growth occurs and is observed.

Note: Other samples of the respiratory tract or related to an associated infectious process may be applied to the present identification system¹⁷⁻²⁰. In the case of nasopharyngeal or pharyngeal exudate samples, the guidelines established for the determination of agents responsible for Respiratory Infections must be followed

The result obtained with this system should not be the sole diagnostic criterion, other clinical criteria or additional tests established by the clinical physician must be taken into account.

6.4 Reading Procedure and Interpretation

Carefully read the present Instructions for use, and carry out the complementary tests indicated.

If, after a 48-hour incubation, no positive reactions are observed in wells 1 to 16, the plate should be incubated and readings should be made every 24 hours, for a total of 5-7 days. All readings are made visually.


7. INTERPRETATION OF RESULTS

The results are made evident by the change in color that occurs in the wells, based on the reactions of the chemical or chromogenic components specifically formulated for each microorganism.

The interpretation is realized through the analysis of the reactions that take place in different wells.

The wells included in the identification panel for the determination of antimicrobial susceptibility have been prepared according to the established lines and have been defined in relation to the antimicrobials suggested by the group of experts consulted in relation to the treatment of the microorganisms included in the the present kit (^{9, 10, 18, 20}). However these can be modified in relation to previous requests or recommendations and in relation to the use of the same. The interpretation of the results is made by analyzing the reactions that occur in various wells (See the Table for the Interpretation of the Results)


Presumptive identification of *Mycoplasma pneumoniae*/*Mycoplasma spp.*




Presence of *Mycoplasma pneumoniae*



Presence of *Mycoplasma pneumoniae*



Presence of *M.pneumoniae*/*Mycoplasma spp.*
(*Mycoplasma penetrans*, *Mycoplasma fermentans*)



Negative results for *Mycoplasma pneumoniae*

ANTIBIOTIC SUSCEPTIBILITY TEST

- A color change to yellow in the field of the antibiogram indicates the presence of bacterial growth and is therefore an index of resistance.
- In the case in which the color of one or more wells turns red, the growth of the test strain will be inhibited and therefore must be interpreted as sensitive.
- Antibiotic concentrations in the wells were prepared according to the guidelines established by CLSI ⁽⁶⁾ for *Mycoplasma pneumoniae* susceptibility testing
- Organisms susceptible to tetracycline are also susceptible to doxycycline ⁽⁶⁾.

8. WARNINGS AND PRECAUTIONS

1. For professional *in vitro* diagnostic use only.
2. Samples should be treated as potentially infectious and tests should only be performed by qualified personnel.
3. Open only the tests necessary to perform the test.
4. Do not use expired devices.
5. Do not use components of any other type of kit.
6. Wear protective clothing and disposable gloves when handling reagents and clinical specimens. Wash your hands thoroughly after performing the test.
7. Do not smoke, drink or eat in areas where specimens or kit reagents are being handled.
8. Dispose of all samples and materials used to perform the test as bio-hazardous waste.
9. Carefully read these Instructions for Use

9. LIMITATIONS

- Samples taken after starting the antimicrobial treatment. A single dose of antibiotic may invalidate the results of the test for detection of agents included in this kit ^{(1) (4)}.
- Inadequate collection and storage of the sample (Paragraph 6.3).
- Inadequate sample (Paragraph 6.3).
- Personnel not employed or trained in microbiology.
- Read these Instructions for Use carefully before carrying out the test in order to avoid errors

10. QUALITY CONTROL

Each batch of MYCO PNEUMO ONE is subjected to a rigorous quality control with different strains. The reference strains are used both for the positive reactions of the different wells as to prove the good functioning of the media formulations of the different wells in case of any non-specific reactions ^{(10) (11) (12)}.

To perform quality control, the following reference strains (including cross-reactions) are recommended:

Mycoplasma buccale ATCC 23636,
Mycoplasma hominis ATCC 23114,
Mycoplasma orale ATCC 23714
Mycoplasma pirum ATCC 25960
Mycoplasma pneumoniae ATCC 15531
Mycoplasma salivarium ATCC 23064
Mycoplasma fermentans ATCC 19989
Haemophilus influenzae ATCC 49247
Streptococcus pyogenes ATCC 19615
Streptococcus agalactiae ATCC 13813
Moraxella catarrhalis ATCC 25238
Staphylococcus aureus ATCC 25923
Enterobacter aerogenes ATCC 13048
Candida albicans ATCC 10231
Pseudomonas aeruginosa ATCC 27853
Klebsiella pneumoniae ATCC 13883
Enterococcus faecalis ATCC 29212

Each laboratory must establish an internal quality control.

11. STORAGE AND CONSERVATION

Store at 2-8 ° C in its original packaging. Do not store near heat sources and avoid extreme temperature variations. In these conditions the product is valid until the expiration date indicated on the label of the primary and secondary container. Do not use after this date. Eliminate if there are signs of deterioration.

12. BIBLIOGRAPHY

1. Atkinson TP, Waites KB. Mycoplasma pneumoniae Infections in Childhood *Pediatr Infect Dis J*. 2014 Jan;33(1):92-4.
2. Atkinson TP. Is asthma an infectious disease? New evidence. [Curr Allergy Asthma Rep](#). 2013 Dec;13(6):702-9.
3. Brunetti L1, Colazzo D, Francavilla R, Tesse R, De Sario V, Lorè M, Armenio L. The role of pulmonary infection in pediatric asthma. [Allergy Asthma Proc](#). 2007 Mar-Apr;28(2):190-3.
4. Chaudhry R1, Ghosh A, Chandolia A. Pathogenesis of *Mycoplasma pneumoniae*: An update. *Indian J Med Microbiol*. 2016 Jan-Mar;34(1):7-16. Bao YX1, Li J1, Tian Y1, Liu QH, Bao J.) Atopy: A risk factor of refractory Mycoplasma pneumoniae pneumonia? Infection-related asthma. *Clin Respir J*. 2015 Dec 19.
5. Chaudhry R1, Ghosh A, Chandolia A. Indian. Pathogenesis of *Mycoplasma pneumoniae*: An update. *J Med Microbiol*. 2016 Jan-Mar;34(1):7-16.
6. Darveaux JI1, Lemanske RF Jr. Infection-related asthma. *J Allergy Clin Immunol Pract*. 2014 Nov-Dec;2(6):658-63.
7. Giavina-Bianchi P, Kalil J. Mycoplasma pneumoniae infection induces asthma onset. [J Allergy Clin Immunol](#). 2016 Apr;137(4):1024-5.
8. Väisänen T, Hämäläinen S, Koivula I, Koskela H. *Mycoplasma pneumoniae* and Stevens-Johnson syndrome. *Respir Res*. 2008 May 30;9:48.
9. Clinical and Laboratory Standard Institute (CLSI), M43-A, Methods for Antimicrobial Susceptibility Testing for Human Mycoplasmas; Approved Guideline, First Edition (2011).
10. Owen B. Spiller. Emerging Pathogenic Respiratory Mycoplasma hominis Infections in Lung Transplant Patients: Time to Reassess its Role as a Pathogen?. *EBioMedicine*. 2017; 19: 8–9.
11. Ng KF et al. Mycoplasma hominis, Ureaplasma urealyticum and Chlamydia trachomatis Respiratory Colonization in Malaysian Preterm Infants. *J Trop Pediatr*. 2017; 2:1092-1093
12. Hagiya H et al. Mycoplasma hominis periaortic abscess following heart-lung transplantation. *Transpl Infect Dis*. 2017;19(3).
13. Patel KK et al. Colonization of paediatric lower respiratory tract with genital Mycoplasma species. *Respirology*. 2011;16(7):1081-7.
14. García C . Isolation of Mycoplasma hominis in critically ill patients with pulmonary infections: clinical and microbiological analysis in an intensive care unit. *Intensive Care Med*. 2007;33(1):143-7.
15. Sánchez, J., *Procedimientos en Microbiología Clínica, Recomendaciones de la Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica*, 2006.
16. Amores J, et al. Viability of *Mycoplasma agalactiae* and *Mycoplasma mycoides* subsp. capri in goat milk samples stored under different conditions. *Vet Microbiol*. 2010;145(3-4):347-350.
17. Boonyayatra S, et al. Effects of storage methods on the recovery of *Mycoplasma* species from milk samples. *Vet Microbiol*. 2010;144(1-2):210-213.
18. Farmacopea Europea. Mycoplasmas. Methods of analysis. 2014; Ed. 8.0, Appendix 2.6.7.
19. Validazione del Sistema A.R.I WELL D-ONE in campioni da I.R.A. Comparazione con metodi tradizionali. A.R.I. WELL D-ONE Technical File, C.P.M. sas, 2013.

COLORIMETRIC TABLE OF RESULTS

COLORIMETRIC REACTION			
INTERPRETATION OF WELL 1			
<p>The well 1 allows the growth of <i>M. hominis</i> and <i>M. pneumoniae /M.spp</i>. The color of the well without inoculation is orange. A red color change in 24/48 hours means growth of <i>M. hominis</i> and, a change in color to yellow after 48 hours means growth of <i>M. pneumoniae / M.spp</i></p>			
1	<i>M.pneumoniae/M.spp</i>	 GIALLO/ARANCIO TRASPARENTE	 ROSSO
		 ROSSO	 GIALLO/ARANCIO TRASPARENTE
INTERPRETATION OF WELLS FROM 2 TO 15			
WELL	DESCRIPTION	POSITIVE	NEGATIVE
2	<i>M.pneumoniae/M.spp</i>	 YELLOW/TRANSPARENT ORANGE	 RED
3	<i>M.pneumoniae/M.spp</i>	 YELLOW/TRANSPARENT GREEN	 RED
4	<i>M.pneumoniae/M.spp</i>	 YELLOW/TRANSPARENT GREEN	 RED
5	<i>M.pneumoniae/M.spp</i>	 YELLOW/TRANSPARENT GREEN	 RED
6	Cross Reaction	 TURBID RED/YELLOW	 TRANSPARENT RED
7	Fungal Growth	 TURBID YELLOW	 SLIGHT YELLOW/TRANSPARENT
8	Bact Growth	 TURBIDITY	 WHITE/SLIGHT YELLOW/ TRANSPARENT
AST			
WELL	DESCRIPTION	RESISTANT	SUSCCEPTIBLE
9	LEV 1	 YELLOW/TRANSPARENT ORANGE	 RED
10	MXF 0,5	 YELLOW/TRANSPARENT ORANGE	 RED
11	TET 2	 YELLOW/TRANSPARENT ORANGE	 RED
12	E 0,5	 YELLOW/TRANSPARENT ORANGE	 RED
13	E 1	 YELLOW/TRANSPARENT ORANGE	 RED
14	AZ 0,5	 YELLOW/TRANSPARENT ORANGE	 RED
15	AZ 1	 YELLOW/TRANSPARENT ORANGE	 RED
INTERPRETATION OF WELL 16			
<p>The well 16 allows the growth of <i>M. hominis</i> and <i>M. pneumoniae /M.spp</i>. The color of the well without inoculation is orange. A red color change in 24/48 hours means growth of <i>M. hominis</i> and, a change in color to yellow after 48 hours means growth of <i>M. pneumoniae / M.spp</i></p>			
16	<i>M. hominis</i>	 RED	 YELLOW/TRANSPARENT ORANGE

<i>M. pneumoniae/M.spp</i>	● ●	YELLOW/TRANSPARENT ORANGE	●	RED
----------------------------	-----	---------------------------	---	-----

Simbols

	Codice	Reference number	Código
	Lotto	Batch number	Número de lote
	Scadenza	Expiration date	Fecha de vencimiento
	Temperatura	Temperature	Temperatura
	Per uso diagnostico <i>in vitro</i>	For <i>in vitro</i> diagnostic use	Para uso de diagnóstico in vitro
	Fragile, manipolare con cura	Fraille, handle with care	Frágil, manejar con cuidado
	Istruzioni per l'uso	Instruction for use	Literatura Interior
	Fabbricante	Manufacturer	Fabricante
	Contenuto sufficiente per <n> test	Content sufficient for <n> tests	Contenido suficiente para <n> pruebas
	Non riutilizzare	Do not re-use	No reutilizar

